

Product No.100225-10- **DEX100**

208.2 LITERS



DIASOL 100225-10- **DEX100**



LIQUID ACID CONCENTRATE FOR BICARBONATE DIALYSIS NON PYROGENIC

Ionic Contribution of Concentrate

Diasol 100225-10- **DEX100**

Diasol liquid acid concentrate is indicated for use in the acute and chronic hemodialysis. It is an accessory to be used with the appropriate hemodialysis machine in a 3 stream mix in the exact prescribed proportion with ANSI/AAMI/ISO 13959:2009 standard water and sodium bicarbonate mix.

SODIUM	100 mEq/L	Final Dialysate Composition: When 1 part acid
POTASSIUM	2.0 mEq/L	concentrate is mixed with 1.72 parts base
CALCIUM	2.5 mEq/L	concentrate (Diasol 10000 Series Bicarbonate)
MAGNESIUM	1.0 mEq/L	and 42.28 parts water, the final ionic
ACETIC ACID	4.0 mEq/L	composition is: sodium 137 mEq/L; chloride
CHLORIDE	105.50 mEq/L	105.50 mEq/L; and net bicarbonate 37 mEq/L.
DEXTROSE	100mg/dL	All other constituents remain the same.

CAUTION: Federal law (USA) Restricts this device to sale by or on order of a physician.

WARNING: For use only with Diasol 10000 Series bicarbonate or equivalent (refer to label). Use of this Concentrate without associated bicarbonate concentrate may cause patient injury or death. **NOT FOR PARENTERAL USE.**

DIASOL, INC.
SAN FERNANDO
CALIFORNIA, 91340
Rev 10/2016

Chemical Composition: Each liter of concentrate contains: 262.98g. Sodium Chloride (NaCl) , 8.28g. Calcium Chloride (CaCl₂·2H₂O) 4.58g. Magnesium Chloride (MgCl₂ ·.6H₂ O), 13.50g Acetic Acid, 6.72g. Potassium Chloride (KCl), 45.00g Dextrose

LOT NO. T03141

MANUFACTURING DATE:
March 14, 2017

NOTE: Verify conductivity and pH prior to use and each time new concentrate is supplied to the machine. Refer to equipment manufacturer's instruction prior to use. Water must meet ANSI/AAMI/ISO standard.

EXPIRATION DATE:
March 14, 2020

RECOMMENDED STORAGE: Avoid excessive heat and protect from freezing.

DO NOT USE IF CONTAINER IS DAMAGED AND/OR SEAL IS BROKEN



B4671002251010055