



CUSTOM DIALYSATE FORMULATIONS!

Diasol/Dryasol

ACID CONCENTRATE

SODIUM BICARB LIQUID AND POWDER

No More Spiking!
No More Lifting and Shaking!
Custom Additives!

- The only “see through drum” which allows you to determine the exact level of concentrate in the drum.
- Large variety of formulations and custom bath such as low/high Ca, high/low K to fit all your patients needs.
- Labels are color coded for easy identification.
- Product code identifies Calcium, Potassium and Magnesium content.
 Ex: 80 2 25-10
 K CA Mg
- Acetate concentrate made upon request.
- Additives

| | | | | |
|----|---------------|---|-----------------------------|----------------------------|
| 1. | 36.83X | — | 8000 Series (Diasol) | 36 Series (Dryasol) |
|----|---------------|---|-----------------------------|----------------------------|

Formulations include: OK -3K, OCA, 2.0CA, 2.5CA, 3.0CA, 3.5CA
 Custom formulation include: 4K, 3.0CA
 All available in Drums, Gallons and Dry Powder
 Sodium Bicarb Powder - Single Treatment, 10 Treatment bags, 1 gallon liquid

| | | | | |
|----|------------|---|-----------------------------|----------------------------|
| 2. | 35X | — | 7800 Series (Diasol) | 35 Series (Dryasol) |
|----|------------|---|-----------------------------|----------------------------|

Formulations include: OK-3K, OCA, 2.0CA, 2.5CA, 3.0CA, 3.5CA
 Custom formulation include: 1.5CA, 1.5K
 All available in Drums, Gallons and Dry Powder
 Sodium Bicarb Powder - Single Treatment, 10 Treatment and 25 Treatment bags, 1 gallon liquid

| | | | | |
|----|------------|---|------------------------------|----------------------------|
| 3. | 45X | — | 10000 Series (Diasol) | 45 Series (Dryasol) |
|----|------------|---|------------------------------|----------------------------|

Formulations include: 0K-3K, 0CA, 2.0CA, 2.5CA, 3.0 CA, 3.5CA
 Custom formulation include: 1.0CA
 All available in Drums, Gallons and Dry Powder
 Sodium Bicarb Powder - Single, 10, 12, 25 Treatment bags, 1 gallon liquid

diasolinc@gmail.com
800/366-0546 • 818/838-7077 • Fax 818/838-7077
www.dialysissupply.com
Now Accepting Visa and Master Card
New! 3 Distribution Centers: San Fernando and Stockton, CA and Phoenix, AZ

Liquid and Dry Concentrate Product Line

| <u>Formulations</u> | <u>Na</u> | <u>K</u> | <u>Ca</u> | <u>Mg</u> |
|--------------------------|-----------|----------|-----------|-----------|
| Liquid/Dry Powder | | | | |
| 100000-75/45000-75 | 101 | 0.0 | 0.0 | .75 |
| 100025-75/45025-75 | 101 | 0.0 | 2.5 | .75 |
| 100110-75/45110-75 | 101 | 1.0 | 1.0 | .75 |
| 100125-75/45125-75 | 101 | 1.0 | 2.5 | .75 |
| 100130-75/45130-75 | 101 | 1.0 | 3.0 | .75 |
| 100210-75/45210-75 | 101 | 2.0 | 1.0 | .75 |
| 100225-75/45225-75 | 101 | 2.0 | 2.5 | .75 |
| 100230-75/45230-75 | 101 | 2.0 | 3.0 | .75 |
| 100325-75/45325-75 | 101 | 3.0 | 2.5 | .75 |
| | | | | |
| 80000-10/36000-10 | 81 | 0.0 | 0.0 | 1.0 |
| 80025-10/36025-10 | 81 | 0.0 | 2.5 | 1.0 |
| 80030-10/36030-10 | 81 | 0.0 | 3.0 | 1.0 |
| 80125-10/36125-10 | 81 | 1.0 | 2.5 | 1.0 |
| 80130-10/36130-10 | 81 | 1.0 | 3.0 | 1.0 |
| 80220-10/36220-10 | 81 | 2.0 | 2.0 | 1.0 |
| 80225-10/36225-10 | 81 | 2.0 | 2.5 | 1.0 |
| 80235-10/36235-10 | 81 | 2.0 | 3.5 | 1.0 |
| 80325-70/36325-70 | 81 | 3.0 | 2.5 | .70 |
| 80330-10/36330-10 | 81 | 3.0 | 3.0 | 1.0 |
| | | | | |
| 78000-10/35000-10 | 105 | 0.0 | 0.0 | 1.0 |
| 78025-10/35025-10 | 105 | 0.0 | 2.5 | 1.0 |
| 78125-10/35125-10 | 105 | 1.0 | 2.5 | 1.0 |
| 78225-10/35225-10 | 105 | 2.0 | 2.5 | 1.0 |
| 78235-10/35235-10 | 105 | 2.0 | 3.5 | 1.0 |
| 78325-10/35325-10 | 105 | 3.0 | 2.5 | 1.0 |

This is the ionic concentration of the acid solution prior to adding Sodium Bicarbonate mix.

If you need a formula that is not listed here — please call, we will make it for you!

800/366-0546 • 818/838-7077 • Fax 818/838-7007

www.dialysissupply.com

Now Accepting Visa and Master Card

New! 3 Distribution Centers: San Fernando and Stockton, CA and Phoenix,